Multiresolution registration of whole IHC slide images to evidence and quantify virtual colocalization of tissue-based biomarkers



DIAPath, Center for Microscopy and Molecular imaging (www.cmmi.be, Gosselies, Belgium)

Whole slide imaging (WSI) using a NanoZoomer 2.0-HT is We evaluated image registration accuracy by measuring the disan essential step of the integrated solution which is developed tances between pairs of control points (located on each image) since 2010 by DIAPath, the digital pathology platform of the which should match if the images were perfectly superimposed. CMMI (Gosselies, Belgium), ensuring the analysis and the Our approach shows accuracy levels compatible with biomarvalidation of tissue-based biomarkers. This solution combines immunohistochemistry (IHC), special staining, WSI, image processing and data analysis. It is based on standardized laboratory procedures and quality controls ensuring reproducibility and traceability. Extracting relevant information from actors involved in complex of registration accuracy into account and validated our combiological processes (such as cancers or treatment responses) re- plete procedure by comparison to colocalization information quires targeting different antigens simultaneously. Multichro- obtained by means of double staining (with different cell locamogenic (brightfield) IHC is used for multiple antigen labeling tions, Fig. 2). on the same slide but suffers from limitations. In addition to antibody cross-reaction problems, this approach also prevents Because registration accuracy was lower for tissue without histhe analysis of proteins expressed in the same cellular compart- tological structures (e.g. the distances between control points ment (because of color merging or masking). We thus deve- being at most 80 µm for serial brain tumor slides), we developed loped an alternative which consists in applying standard IHC a sequential IHC technique applied on the same slide (Fig. 1, on consecutive (also labeled "serial") tissues sections to target right-hand panel). This "Sequential Immunoperoxidase Labeling different proteins (see Fig. 1, left-hand panel). These slides were and Erasing" (SIMPLE) method is based on cycles of staining/ then digitalized and superimposed by a registration algorithm. digitization/erasing, where after IHC staining and slide digiti-Our approach uses the fact that histological structures (such as zation, staining is "washed" through an antibody elution techglands or epithelium) are often well conserved across a few se- nique. We improved the original SIMPLE method by including rial slides because the slide thickness obtained with microtomy a new elution methodology which preserves tissue and antigen is 3 to 5 μ m (i.e. smaller than the mean size of human cells).

ker colocalization characterization. Indeed, on serial slides the distances between control points were evaluated to be at most 20 µm in presence of histological structures in the tissue (e.g. in colonic tissue, see Fig. 2). We then implemented a method to extract biomarker colocalization measurements taking the level

2) Registration of pairs of core images

3) Colocalization map overlap



To superimpose the serial slide images we developed an effition images (showing the whole slides to register in small-sized points [1]. images such as 4,000 by 3,000 pixels). Briefly, this step enables

epitopes for the next staining step. We successfully used this approach to identify antigens expressed in the same cellular compartment of high-grade glioma samples. We tested our recient 2-step registration algorithm [1]. The first, low-resolution, gistration method on the virtual slides so obtained and achieved registration step is applied on the 1X equivalent magnifica- very good results, i.e. about 5 µm of distance between control to roughly superimpose two virtual slides on the basis of the slice Afterwards we adapted our registration process to tissue miarray fashion. As expected, we obtained very accurate registration results on serial TMA slides, even for unstructured tissue (i.e. distance between control points less than $10 \,\mu\text{m}$) [2]. Fig. 3 shows the different steps of the process requiring (1) the correct identification of the circular samples (from the cores) on each TMA image, (2) the registration of core images extracted from serial slide images and (3) the computation of a colocalization map on small squared tiles from which a global colocalization measurement (adapted from immunofluorescence) can be computed [1,2]. A first (not illustrated) step consists in computing in each tile the value of a staining feature characterizing each IHC (1) DIAPath, Center for Microscopy and Molecular Imaging, marker, e.g. the labeling index (LI), which is the percentage of positive (i.e. DAB-stained) tissue area [3]. The colocalization map evidences local staining overlap, which increases when both LI increases in a tile (see Fig. 3.3). The value of the global values between 0 and 1.

Figure 3: IHC staining colocalization on serial TMA slides showing IGF-I and IGFBP2 IHC expression in colonic tissue. After 1) core identification on each TMA image, 2) core images are extracted and registered to compute 3) staining overlap map and index (r0). In the map blue indicates no overlap whereas yellow to red show increasing overlap levels.

References:

[1] X. Moles Lopez, et al. Registration of whole immunohistochemical slide images: an efficient way to characterize biomarker colocalization. J Am Med Inform Assoc 22(1):86–99, 2015. [2] Y.-R. Van Eycke, et al. High-throughput analysis of tissue-based biomarkers in digital pathology. In: 2015 37th Annu. Int. Conf. IEEE Eng. Med. Biol. Soc. IEEE, Aug. 2015, pp. 7732-7735.

[3] C. Decaestecker, et al. Requirements for the valid quantification of immunostains on tissue microarray materials using image analysis. Proteomics 9(19):4478–94, 2009.

outlines, holes and/or large histological structures. After this first croarray (TMA) slides for high-throughput analysis of biomarstep a series of high-resolution registrations were independent- ker colocalization. Registration is easier for these slide images ly performed on the set of high-resolution (i.e. 20X equivalent) because TMA assembles cores extracted from tissue blocks in fields of view constituting a whole slide image.



Figure 1: IHC on serial tissue slides (left-hand panel) and overlap index (r0) is the sum of the tile contributions and takes SIMPLE technique applied on the same slide (right-hand panel). Both provide slide images to superimpose for colocalizing different antigen expressions.



1) Core identification on a virtual TMA slide



Yves-Rémi Van Eycke (1,2), Justine Allard (1), Cédric Balsat (1), Isabelle Salmon (1,3), Christine Decaestecker (1,2)

Université Libre de Bruxelles (ULB), Gosselies, Belgium

(2) Laboratories of Image, Signal processing & Acoustics, ULB, Brussels, Belgium

(3) Department of Pathology, Erasme Hospital, ULB, Brussels, Belgium

For further information, please contact: amarcowycz@hamamatsu.be